

Quality. Simplicity. Control.

MediCal QAWeb is the industry's first, and most deployed, web-based automated quality assurance solution. Used in more than 3,500 facilities worldwide, the all-inclusive, secured system is your guarantee for consistent image quality, compliance and uptime of all PACS display systems throughout your facility.

Take full control of your hospital's assets by installing QAWeb today.





Do more
with less

Because their quality and uptime are critical for productivity and accuracy, diagnostic systems need to be perfectly reliable. Hospitals are also looking to find creative financial solutions to reduce costs while still optimizing their resources. How do you respond to these seemingly disparate challenges?

MediCal QAWeb is the answer.

Seamless Compatibility

MediCal QA Web is fast and easy to implement, requires no additional hardware and works with all your displays—**Barco and non-Barco alike**

Ultimate Peace of Mind

Calibration and QA tasks are performed automatically without interrupting the viewing process—ensuring 24/7 uptime

Effortless Asset Management

With data analytics you get a comprehensive overview of your workstations and can proactively plan, justify and prioritize capital expenses for hardware

MediCal QAWeb enables us to sustain productivity levels despite fewer staff, while increasing operational efficiency across the hospital.

Cathy Taylor, Manager of the Clinical Applications Team, Gwinnett Medical Center

An all-inclusive and secure solution

Remote and local access

Through inclusive enterprise management users can easily log in from anywhere and instantly check the status of every display or undertake all necessary corrective actions if an issue arises

Security at every level


The system—available via a secured connection—makes use of a central and HIPAA compliant server that provides users with access to all relevant QA data from anywhere at any time

Straightforward compliance

Meeting the changing compliance standards of internationally accepted medical guidelines is made easy with automated reporting and automated DICOM calibration



Enterprise Solution: Multiple Hospitals | Multiple Locations



MediCal QAWeb Software Licenses

Feature / Software License		Essential	Clinical	Premium
		Available out of the box, free of charge	Additional level of support with full reporting features	Extended central & remote QA and asset management functionality with full reporting
(Pre)Software update notifications		Yes	Yes	Yes
Software upgrades			Yes	Yes
Access to knowledge base and issue escalation		Yes	Yes	Yes
Remote support		Yes	Yes	Extended
QA Profile		Essential	Clinical	Premium
Asset Management	Central hardware fleet management	Basic	Extended	Extended
	Central hardware & software asset management		Yes	Yes
	Issue notification & knowledge-based solutions		Yes	Yes
QA Management	Central calibration control	Yes*	Yes*	Yes
	Central scheduling control	Yes	Yes	Extended
	QA by functional hardware & software checks	Basic	Basic	Extended
	QA issue notification & knowledge-based solutions		Yes	Yes
	Visual Confidence Test		Yes	Yes
Report Management	General reports		Yes	Yes
	Activity reports	Yes	Yes	Yes
	Budget reports		Yes	Yes
	JCAHO reports		Yes	Yes
	History reports		Yes	Yes

All service levels are available after registration to the online MediCal QAWeb Server

* Remote calibration is initiated after applying a new policy